



Ultrafast – Simply Brilliant



InstantBlue™ Ultrafast protein stain



Ultra Quick

Results within 15 minutes

Easy

One step protocol

- No washing
- No fixing
- No heating
- No de-staining

Sensitive


≥5ng protein per band

Enabling


Mass Spec compatible
Edman Sequencing

Safe

Non toxic
No Methanol
Sink disposal



Innovators of Protein Technologies



Does your gel stain tick all the boxes?

	InstantBlue™ (Expedeon)	SimplyBlue™ (Invitrogen)	GelCode™ Blue (Pierce)	Home-made (Anyone)
No need to wash before staining	✓	3 x 5 X minutes	3 x 5 X minutes	Typically 3 x 5 X minutes
No need to fix before staining	✓	✓	X ³	✓
No need to heat/microwave	✓	X ¹	✓ ⁴	✓ ¹
No need to de-stain	✓	X	X	X
No risk of over-staining if left too long	✓	X ²	X	X
Non toxic and easily disposed of	✓	✓	X	X
True single step staining	✓	X	X	X

¹ Recommended, otherwise requires 1 or more hours to stain ² Sensitivity decreases after 3 hours staining ³ If gels run with MOPS or MES buffer systems ⁴ Suggested to provide faster, 30 minute staining but sensitivity reduced



InstantBlue™ Ultrafast protein stain

For further information go to
www.expedeon.com

Expedeon Ltd
Unit 1A, Button End, Harston,
Cambs., CB22 7GX United Kingdom
t: +44 (0)1223 873364 f: +44 (0)1223 873371

 **Westburg** | (+31) 33 494 6666
Belgium 0800 - 19815
Head Office: P.O.Box 214 | 3830 AE Leusden | The Netherlands
fax (+31) 33 495 1222 | info@westburg.eu | www.westburg.eu